



Emergencies Happen!

Dynarex is there for you!

Thousands of EMS products for all emergency situations.



Add These New EMS Products to Your Offering!

Heavy Duty Black Body Bags



Post-Mortem Bag Kits



Cardboard Splints with Foam



Also Good for First Aid/Kit Packers!

Mini Scissors



Wrapped in an individual poly bag.



Bandage Shears



Wrapped in an individual poly bag.



See Reverse Side for More Product Details



The Brand To Know for EMS Products

Heavy Duty Black Body Bag 4776



Reorder No.	Description	Size	Quantity
4776	Heavy Duty Black Body Bag	36"x 90"	10/Cs

Dynarex heavy duty black body bag comes with six handles for easy transport and has an envelope style closure. Includes three toe tags and hazard labels. 12 mil thickness. Double zipper construction for easy closure.

Post-Mortem Bag Kits 4772-4775



Reorder No.	Description	Size	Quantity
4772	Adult Body Bag Kits	36"x 90"	10/Cs
4773	Adult Body Bag Kits	48"x 90"	10/Cs
4774	Pediatric Body Bag Kits	28"x 48"	10/Cs
4775	Bariatric Body Bag Kits	72"x 96"	5/Cs

Made of lighter material, this body bag kit comes in four sizes, has an envelope style closure and includes three toe tags and hazard labels. 7 mil thickness. Double zipper construction for easy closure.

Cardboard Splint with Foam 3624-3626



Reorder No.	Description	Size	Quantity
3624	Cardboard Splint	9"x12"	30/Cs
3625	Cardboard Splint	9"x18"	30/Cs
3626	Cardboard Splint	14"x24"	30/Cs

Lightweight disposable splints with 3/8" foam. Used in EMS settings for stabilizing fractures of extremities. Conforms easily to limbs. Holes can be used to secure splints to injured extremity. This product is made from cardboard, foam and glue. It does not contain metal and is radiolucent.

Also Good for First Aid/Kit Packers!

Mini Scissors 4190



Reorder No.	Description	Size	Quantity
4190	Mini Scissors	3.5"	5/100/Cs

Stainless steel construction with polyethylene handles. Perfect for first aid kits and emergency situations. Can cut gauze, tape and other small items. **Wrapped in an individual poly bag.**

Bandage Shears 4180-4187



Reorder No.	Description	Size/Color	Quantity
4180	Bandage Shears	Black	4/50/Cs
4181	Bandage Shears	Blue	4/50/Cs
4182	Bandage Shears	Magenta/Pink	4/50/Cs
4183	Bandage Shears	Orange	4/50/Cs
4184	Bandage Shears	Purple	4/50/Cs
4185	Bandage Shears	Red	4/50/Cs
4186	Bandage Shears	Teal	4/50/Cs
4187	Bandage Shears	Yellow	4/50/Cs

Bandage shears come in multiple colors. Stainless steel construction with one serrated edge. Strong enough to cut clothing and seat belts. **Wrapped in an individual poly bag.**

MANY MORE DYNAREX EMS PRODUCTS TO CHOOSE FROM



Visit www.dynarex.com

Find us on:

