

## Handy Sandy Sand Bags From Morrison Medical

- Tough nylon reinforced vinyl exterior.
- Dielectrically sealed to prevent leakage or escape of dust.
- Sewn in weight labels and color coded by weight for easy identification.
- Complete with convenient handle.

### Handy Sandy Sand Bags

*Please Call for Pricing*



*Ask us about screen printing for additional identification, as "OR" is shown above.*

**Item #0360:**

3 lbs. Handy Sandy Sand Bag, Yellow, 4½" x 9½", each

**Item #0360-4:**

3 lbs. Handy Sandy Sand Bag, Yellow, 4½" x 9½", case of 4

**Item #0370:**

5 lbs. Handy Sandy Sand Bag, Orange, 5½" x 12", each

**Item #0370-4:**

5 lbs. Handy Sandy Sand Bag, Orange, 5½" x 12", case of 4

**Item #0380:**

7 lbs. Handy Sandy Sand Bag, Red, 6½" x 15½", each

**Item #0380-4:**

7 lbs. Handy Sandy Sand Bag, Red, 6½" x 15½", case of 4

**Item #0390:**

10 lbs. Handy Sandy Sand Bag, Blue, 6½" x 17", each

**Item #0390-4:**

10 lbs. Handy Sandy Sand Bag, Blue, 6½" x 17", case of 4

*For more information, call 800-438-6677,  
or visit [www.MorrisonMed.com](http://www.MorrisonMed.com).*

